

OROBOROS O2k-Manual

Mitochondrial Physiology Network 21.16(02):1-8 (2016)
Version 02: 2016-08-15 ©2016 OROBOROS
Updates: http://wiki.oroboros.at/index.php/MiPNet21.16_DatLab_7



DatLab 7: innovations

Gradl L¹, Gnaiger E^{2,3}, Capek O³, Plattner C⁴



OROBOROS INSTRUMENTS

Schöpfstr 18, A-6020 Innsbruck, Austria
Email: instruments@oroboros.at
www.oroboros.at


Overview




DatLab 7 is a next step in high-resolution respirometry (HRR) incorporating user-friendly features on quality control, documentation, and traceability of measurement with the O2k (MitoFit). Conceptual unification is obtained in calibrations and flux analysis of O2k-MultiSensor channels, including normalization and baseline correction. Graphs are further improved for real-time display and publication.

1. General

1.1. User




For better documentation and traceability, a user name has to be entered or selected  after starting DatLab. The current user is displayed in the O2k signal line. User names are connected with personal graph layouts (see [Layout](#) menu).

Change user: When DatLab is connected to the O2k, the user can be changed under [File](#) \ [Change user](#).

Manage users: Existing users can be renamed and deleted (with all connected graph layouts) under [File](#) \ [Manage users](#).

1.2. Menus



New, user-friendly structure of menus: **Flux/Slope** as a separate window with new features (see Section 6).

1.3. Hyperlinks in DatLab windows

In many DatLab windows, a link to the OROBOROS website leads directly to context-sensitive help.

1.4. Event names

Event names are shown on top of the graph to avoid text overlaps with mark names.

1.5. Illumination

“Light” has been renamed to “Illumination” to distinguish it from the LED light introduced into the chambers in the O2k-Fluorometer.

2. Oxygraph-2k

In the Oxygraph-2k menu, the “Show channel”-functions have been moved and an additional submenu [O2k configuration] added for a better overview.

Oxygraph-2k	TIP2k	Experiment
O2k control		F7
O2k configuration		
Stirrer A on/off		F11
Stirrer B on/off		F12
Stirrer test		F9
Illumination on/off		F10
Manage setups		

2.1. Data recording interval The data recording interval can no longer be changed when the O2k is running. It has to be set when the O2k control window pops up at the beginning of an experiment.

2.2. O2k configuration window

 A screenshot of the 'O2k configuration' dialog box. It has a title bar 'O2k configuration'. Inside, there are several sections:

- 'O2k serial number' with a text input field.
- 'Power-O2k' with a 'P' label and a text input field.
- 'Chamber' with 'A' and 'B' labels and two text input fields.
- 'Oxygen, O2' section with 'Oxygen sensor #' and 'Channel label' (containing 'O2') for both chambers.
- 'Amperometric, Amp' section with a checked checkbox, 'Amp sensor #' and 'Channel label' (containing 'Amp') for both chambers.
- 'Potentiometric, pX' section with a checked checkbox, 'pX Electrode #' and 'pX Reference electrode #' text input fields, and 'Channel label' (containing 'pX') for both chambers.
- At the bottom, there is a checked checkbox 'Skip configuration at reconnect', a 'Cancel' button, and an 'OK' button. A link 'MitoPedia: O2k configuration' is also present.

In **Oxygraph-2k\O2k configuration**, channels (Amperometric, Amp , Potentiometric, pX) can be switched on/off by selecting the according tick box.


Power-O2k (P#) and numbers of Oxygen sensors, Amp sensors, pX electrodes and pX reference electrodes can be entered and edited. They are displayed in the **O2k control** window.

O2k configuration pops up automatically after connecting to the O2k.

- Skip configuration at reconnect If selected, **O2k control** pops up immediately when re-connecting to the O2k.


3. TIP2k

3.1. Insert/Append



A new button **Append** improves user-friendly programming of the TIP2k, for adding a new line at the end of the programme. **Insert** adds a new line at the selected position.

3.2. TIP2k configuration



The tab **Configuration** accommodates the TIP2k series and Microsyringe type. Based on this information, DatLab checks the total volume programmed in the **Control** tab and warns the user if the approximate maximum operation volume is exceeded.

4. Experiment

4.1. Edit experiment

The Edit experiment window is complemented with new entry boxes:

File recorded by (*automatically filled*)

Protocol: Protocol name

Sample type: former "Sample"

Cohort

Sample code

Sample number

Subsample number

Protocol and sample number are shown next to the graph (right side) above the displayed channel names.

⌘ A double-click on Protocol name or Sample number opens the **Edit experiment** window.




4.2. Experimental log

Experimental log replaces the name **View protocol**.


5. Calibration

5.1. Calibration of pX channel

- The new **pX calibration** window allows multipoint calibration, comparable to the **Amp calibration** window.
- Units for the calibrated pX signal and pX slope were updated, eg pH and mpH.



6. Flux/Slope



The previous tab **Slope** in the **Calibration** window has been moved to the new window **Slope configuration** selected for each O2k chamber and channel under the **Flux/Slope** menu. New features are added and selecting with the according radio button or tick box.

- ⊙ **Flux per volume**
- ⊙ **Flux per mass** or **Flow per cells**, depending on entries in the **Edit experiment** window, as displayed.
- ⊙ **Flux control ratio, FCR**
- ☑ **Baseline correction** (new)
- ☑ **Background correction**
- ▼ **Slope smoothing**
- ⊙ **Stoichiometry** (new): default -1 for O2 and +1 for Amp and pX.

7. Graph

7.1. Select plots – Graph layout

The **Graph \ Select plots** window has been restructured and renamed as **Graph layout**.

In the tab **Plots**, channels are separated in different tabs. For each **Graph**▼, plots are selected with tick boxes for the Y-axes Y1 or Y2.

The axis default labels are changed automatically, according to the channel and unit. Signal units selected in **Plots** are independent of the units used for calibration.

To display the **Raw signal**, the tick box on the right side of the according channel has to be selected.

7.2. Custom label

In DatLab 7, users have the possibility to change the axis labels in [Graph/Select plots]. Two lines are available for entering a name and the according unit.

8. Layout

8.1. Layout structure

A new hierarchy for DatLab graph layouts is implemented by separation of the following **Layout categories**:


Reference layouts are provided with the installation package of DatLab 7. They cannot be changed or deleted, but it is possible to edit a reference layout and save it under a new name.

All users These layouts can be used and edited by all users.

User: Name These personal layouts can only be accessed and edited by the User defined upon starting DatLab.

For every layout category, the following **Layout types** are provided:

- ▶ **O2k-Core**
- ▶ **O2 & Amp**
- ▶ **O2 & pX**
- ▶ **Other**




8.2. Most recently used layouts

The last five layouts used are displayed in the menu **Layout**.

9. Marks

9.1. Mark statistics – Copy to clipboard options



Experimental details are copied to clipboard as default in addition to Mark statistics information. Deselect the tick box for using older versions of DatLab-Excel templates.

Traceability in DatLab-Excel templates is the default to show and export Flux/Slope values as uncorrected slopes (positive or negative) for traceability, where corrections and normalization are transparent on the basis of **Experimental details**. Deselect the tick box to show and copy the current values (e.g. flux per mass).

9.2. Mark statistics – Select

- Median is exported as default (*new*) instead of the averages of values in the marks.
- Range is an additional new export option.

Select

Median

Average

Range

Maximum

Minimum

9.3. Edit mark information window

The **Edit mark information** window (opened by a click on the mark bar) is complemented with three new functions:

- Delete points** deletes all data points in the plot within the selected mark.
- Interpolate points** interpolates all data points in the plot within the selected mark section from the last previous to the first following data point.
- Restore points** reverses data deletion or interpolation in a signal plot within the marked section.
- Recalc slope** recalculates the slope for an entire Flux/Slope plot, thus reversing all data deletions or interpolations in this plot.

Edit mark information

Start 00:27:22

Stop 00:28:05

N Points 44

Average 852.5082

Name

Value

Comment

[MitoPedia: Marks](#)

Delete points

Interpolate points

Recalc. slope

Cancel OK



The joy of success is the next step

Supplement

A. Author contributions

- ¹software security networks – ssn, Innsbruck. LG is the programmer of DatLab 7 and all previous DatLab versions and contributed to the concept of quality control in DatLab 7. External consultant of the MitoFit project.
- ²Universitätsklinik für Visceral-, Transplantations- und Thoraxchirurgie, D. Swarovski Forschungslabor, DSL, Medizinische Universität Innsbruck; ³OROBOROS INSTRUMENTS. EG contributed to the concept of DatLab 7 and all previous DatLab versions, edited the final version of this MiPNet publication, and edited complementary websites (help) on the OROBOROS-Bioblast homepage. MitoFit project leader.
- ³OROBOROS INSTRUMENTS. OC was mainly responsible for trouble shooting of test versions from DatLab 6 to DatLab 7 and contributed to writing this MiPNet publication.
- ⁴Biozentrum Innsbruck - Sektion für Bioinformatik, SBI, Medizinische Universität Innsbruck. CP wrote this MiPNet publication and edited complementary websites (help) on the OROBOROS-Bioblast homepage (MitoFit project partner Univ.-Prof. Dr. Zlatko Trajanoski).

B. Acknowledgements



Contribution to K-regio project MitoFit. The project MitoFit is funded by the Land Tirol within the programme K-Regio of Standortagentur Tirol. www.mitofit.org

We acknowledge contributions by all OROBOROS-team members, particularly Carolina Doerrier, Verena Laner, Elisabeth Hiller, Stephanie Driescher, Valentina Dikova, and Zuzana Sumbalova.